



# DDX3Y Polyclonal Antibody

<b>Catalog No</b>	YP-Ab-01643
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse
<b>Applications</b>	IHC;IF;ELISA
<b>Gene Name</b>	DDX3Y
<b>Protein Name</b>	ATP-dependent RNA helicase DDX3Y
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human DDX3Y. AA range:41-90
<b>Specificity</b>	DDX3Y Polyclonal Antibody detects endogenous levels of DDX3Y protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	IHC: 1/100 - 1/300. ELISA: 1/20000.. IF 1:50-200
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	DDX3Y; DBY; ATP-dependent RNA helicase DDX3Y; DEAD box protein 3; Y-chromosomal
<b>Observed Band</b>	
<b>Cell Pathway</b>	Cytoplasm . Nucleus . Shuttles between the nucleus and the cytoplasm in an XPO1-dependent manner.
<b>Tissue Specificity</b>	Widely expressed at the mRNA level, with highest levels in testis (PubMed:9381176). Testis-specific (at protein level). Expressed predominantly in spermatogonia, but occasionally detected in some pre-leptotene/leptotene spermatocytes (PubMed:15294876).
<b>Function</b>	disease:DDX3Y is deleted in severe non-obstructive hypospermatogenesis [MIM:400042].,function:Probable ATP-dependent RNA helicase. May play a role in spermatogenesis.,similarity:Belongs to the DEAD box helicase family. DDX3/DED1 subfamily.,similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,subcellular location:Shuttles between the nucleus and the cytoplasm in an XPO1-dependent manner.,subunit:May interact with TDRD3.,tissue specificity:Testis-specific. Expressed predominantly in spermatogonia.,
<b>Background</b>	The protein encoded by this gene is a member of the DEAD-box RNA helicase family, characterized by nine conserved motifs, included the conserved Asp-Glu-Ala-Asp (DEAD) motif. These motifs are thought to be involved in ATP



binding, hydrolysis, RNA binding, and in the formation of intramolecular interactions. This protein shares high similarity to DDX3X, on the X chromosome, but a deletion of this gene is not complemented by DDX3X. Mutations in this gene result in male infertility, a reduction in germ cell numbers, and can result in Sertoli-cell only syndrome. Pseudogenes sharing similarity to both this gene and the DDX3X paralog are found on chromosome 4 and the X chromosome. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Oct 2014],


**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using DDX3Y Antibody. The picture on the right is blocked with the synthesized peptide.