





AKR1CL1 Polyclonal Antibody

Catalog No	YP-Ab-02487
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IF;ELISA
Gene Name	AKR1CL1
Protein Name	Aldo-keto reductase family 1 member C-like protein 1
Immunogen	The antiserum was produced against synthesized peptide derived from human AKR1CL1. AA range:80-129
Specificity	AKR1CL1 Polyclonal Antibody detects endogenous levels of AKR1CL1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AKR1CL1; Aldo-keto reductase family 1 member C-like protein 1
Observed Band	37kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Testis,
Function	similarity:Belongs to the aldo/keto reductase family.,
Background	similarity:Belongs to the aldo/keto reductase family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunofluorescence analysis of A549 cells, using AKR1CL1 Antibody. The picture on the right is blocked with the synthesized peptide.