




ARRDC2 Polyclonal Antibody


Catalog No	YP-Ab-03719
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;IHC;IF;ELISA
Gene Name	ARRDC2
Protein Name	Arrestin domain-containing protein 2
Immunogen	The antiserum was produced against synthesized peptide derived from human ARRDC2. AA range:141-190
Specificity	ARRDC2 Polyclonal Antibody detects endogenous levels of ARRDC2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARRDC2; PP2703; Arrestin domain-containing protein 2
Observed Band	44kD
Cell Pathway	plasma membrane,cytoplasmic vesicle,
Tissue Specificity	Brain,Synovial cell,
Function	similarity:Belongs to the arrestin family.,
Background	similarity:Belongs to the arrestin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images



Western Blot analysis of various cells using ARRDC2 Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4° overnight). High-pressure and temperature Tris-EDTA, pH 8.0 was used for antigen retrieval. Negative control (right) obtained from antibody was pre-absorbed by immunogen peptide.



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using ARRDC2 Antibody. The picture on the right is blocked with the synthesized peptide.