



SH3RF3 Polyclonal Antibody

Catalog No	YP-Ab-04207
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	SH3RF3
Protein Name	SH3 domain-containing RING finger protein 3
Immunogen	Synthesized peptide derived from the C-terminal region of human SH3RF3.
Specificity	SH3RF3 Polyclonal Antibody detects endogenous levels of SH3RF3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SH3RF3; POSH2; SH3MD4; SH3 domain-containing RING finger protein 3; Plenty of SH3s 2; SH3 multiple domains protein 4
Observed Band	92kD
Cell Pathway	
Tissue Specificity	Brain,Spleen,
Function	domain:The third SH3 domain mediates interaction with PAK2.,similarity:Belongs to the SH3RF family.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 4 SH3 domains.,subunit:Interacts with PAK2.,
Background	domain:The third SH3 domain mediates interaction with PAK2.,similarity:Belongs to the SH3RF family.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 4 SH3 domains.,subunit:Interacts with PAK2.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images