



# Collagen XV α1 (Cleaved-Asn1212) rabbit pAb

|                                  |                                                                                                                                                     |
|----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                | YP-Ab-04380                                                                                                                                         |
| <b>Isotype</b>                   | IgG                                                                                                                                                 |
| <b>Reactivity</b>                | Human;Mouse                                                                                                                                         |
| <b>Applications</b>              | WB; ELISA                                                                                                                                           |
| <b>Gene Name</b>                 | COL15A1                                                                                                                                             |
| <b>Protein Name</b>              | Collagen XV α1 (Cleaved-Asn1212)                                                                                                                    |
| <b>Immunogen</b>                 | Synthesized peptide derived from human Collagen XV α1 (Cleaved-Asn1212)                                                                             |
| <b>Specificity</b>               | This antibody detects endogenous levels of Human,Mouse Collagen XV α1 (Cleaved-Asn1212, protein was cleaved amino acid sequence between 1211-1212 ) |
| <b>Formulation</b>               | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                             |
| <b>Source</b>                    | Polyclonal, Rabbit,IgG                                                                                                                              |
| <b>Purification</b>              | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                           |
| <b>Dilution</b>                  | WB 1:1000-2000 ELISA 1:5000-20000                                                                                                                   |
| <b>Concentration</b>             | 1 mg/ml                                                                                                                                             |
| <b>Purity</b>                    | ≥90%                                                                                                                                                |
| <b>Storage Stability</b>         | -20°C/1 year                                                                                                                                        |
| <b>Synonyms</b>                  | Collagen alpha-1(XV) chain [Cleaved into: Endostatin (Endostatin-XV;Related to endostatin;Restin)]                                                  |
| <b>Observed Band</b>             | 20 150kD                                                                                                                                            |
| <b>Cell Pathway</b>              | Secreted, extracellular space, extracellular matrix .                                                                                               |
| <b>Tissue Specificity</b>        | Expressed predominantly in internal organs such as adrenal gland, pancreas and kidney.                                                              |
| <b>Function</b>                  |                                                                                                                                                     |
| <b>Background</b>                |                                                                                                                                                     |
| <b>matters needing attention</b> | Avoid repeated freezing and thawing!                                                                                                                |
| <b>Usage suggestions</b>         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                   |



## Products Images