



Eosinophil peroxidase (heavy chain, Cleaved-Val251) rabbit pAb

Catalog No	YP-Ab-04382
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	EPX EPER EPO EPP
Protein Name	Eosinophil peroxidase (heavy chain, Cleaved-Val251)
Immunogen	Synthesized peptide derived from human Eosinophil peroxidase (heavy chain, Cleaved-Val251)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat Eosinophil peroxidase (heavy chain, Cleaved-Val251, protein was cleaved amino acid sequence between 250-251)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Eosinophil peroxidase (EPO;EC 1.11.1.7) [Cleaved into: Eosinophil peroxidase light chain; Eosinophil peroxidase heavy chain]
Observed Band	50 75kD
Cell Pathway	Cytoplasmic granule. Cytoplasmic granules of eosinophils.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images