



NOTCH3 (Cleaved-Asp1628) rabbit pAb

Catalog No	YP-Ab-04396
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	NOTCH3
Protein Name	NOTCH3 (Cleaved-Asp1628)
Immunogen	Synthesized peptide derived from human NOTCH3 (Cleaved-Asp1628)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat NOTCH3 (Cleaved-Asp1628, protein was cleaved amino acid sequence between 1628-1629)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Neurogenic locus notch homolog protein 3 (Notch 3) [Cleaved into: Notch 3 extracellular truncation; Notch 3 intracellular domain]
Observed Band	78 260kD
Cell Pathway	Cell membrane ; Single-pass type I membrane protein.; [Notch 3 intracellular domain]: Nucleus. Following proteolytical processing NICD is translocated to the nucleus.
Tissue Specificity	Ubiquitously expressed in fetal and adult tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images