



PXDC2 Polyclonal Antibody

Catalog No	YP-Ab-05001
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	PLXDC2 TEM7R UNQ2514/PRO6003
Protein Name	Plexin domain-containing protein 2 (Tumor endothelial marker 7-related protein)
Immunogen	Synthesized peptide derived from human protein . at AA range: 50-130
Specificity	PXDC2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	58kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Expressed in the endothelial cells of the stroma but not in the endothelial cells of normal colonic tissue.
Function	function:May play a role in tumor angiogenesis.,similarity:Belongs to the plexin family.,similarity:Contains 1 PSI domain.,subunit:Interacts with CTTN.,tissue specificity:Expressed in the endothelial cells of the stroma but not in the endothelial cells of normal colonic tissue.,
Background	function:May play a role in tumor angiogenesis.,similarity:Belongs to the plexin family.,similarity:Contains 1 PSI domain.,subunit:Interacts with CTTN.,tissue specificity:Expressed in the endothelial cells of the stroma but not in the endothelial cells of normal colonic tissue.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images