



FAKD5 Polyclonal Antibody

Catalog No	YP-Ab-05005
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	FASTKD5 KIAA1792
Protein Name	FAST kinase domain-containing protein 5
Immunogen	Synthesized peptide derived from human protein . at AA range: 500-580
Specificity	FAKD5 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	84kD
Cell Pathway	Mitochondrion matrix, mitochondrion nucleoid . Localizes to mitochondrial RNA granules found in close proximity to the mitochondrial nucleoids. .
Tissue Specificity	Expression detected in spleen, thymus, testis, ovary, colon, heart, smooth muscle, kidney, brain, lung, liver and white adipose tissue with highest expression in heart, smooth muscle, liver and thyroid.
Function	similarity:Belongs to the FAST kinase family.,similarity:Contains 1 RAP domain.,
Background	similarity:Belongs to the FAST kinase family.,similarity:Contains 1 RAP domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

