



ALPK3 Polyclonal Antibody

Catalog No	YP-Ab-05016
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	ALPK3 KIAA1330 MAK
Protein Name	Alpha-protein kinase 3 (EC 2.7.11.-) (Muscle alpha-protein kinase)
Immunogen	Synthesized peptide derived from human protein . at AA range: 730-810
Specificity	ALPK3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	209kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,
Function	function:Kinases that recognize phosphorylation sites in which the surrounding peptides have an alpha-helical conformation. Plays a role in cardiomyocyte differentiation.,similarity:Belongs to the protein kinase superfamily. Alpha-type protein kinase family. ALPK subfamily.,similarity:Contains 1 alpha-type protein kinase domain.,similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,
Background	function:Kinases that recognize phosphorylation sites in which the surrounding peptides have an alpha-helical conformation. Plays a role in cardiomyocyte differentiation.,similarity:Belongs to the protein kinase superfamily. Alpha-type protein kinase family. ALPK subfamily.,similarity:Contains 1 alpha-type protein kinase domain.,similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images