



BMP2K Polyclonal Antibody

Catalog No	YP-Ab-05026
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	BMP2K BIKE HRIHFB2017
Protein Name	BMP-2-inducible protein kinase (BIKe) (EC 2.7.11.1)
Immunogen	Synthesized peptide derived from human protein . at AA range: 740-820
Specificity	BMP2K Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	127kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Epithelium,Ovary,Synovial membrane tissue,Testis,Thyroid,
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,function:May be involved in osteoblast differentiation.,PTM:Autophosphorylated.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family.,similarity:Contains 1 protein kinase domain.,
Background	This gene is the human homolog of mouse BMP-2-inducible kinase. Bone morphogenic proteins (BMPs) play a key role in skeletal development and patterning. Expression of the mouse gene is increased during BMP-2 induced differentiation and the gene product is a putative serine/threonine protein kinase containing a nuclear localization signal. Therefore, the protein encoded by this human homolog is thought to be a protein kinase with a putative regulatory role in attenuating the program of osteoblast differentiation. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images