



ITGA9 Polyclonal Antibody

Catalog No	YP-Ab-05110
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	ITGA9
Protein Name	Integrin alpha-9 (Integrin alpha-RLC)
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	ITA9 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	113kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	In airway epithelium, in the basal layer of squamous epithelium, and in smooth muscle, skeletal muscle, and hepatocytes. Abundantly expressed in fetal lung and lung cancers.
Function	function: Integrin alpha-9/beta-1 is a receptor for VCAM1, cytotactin and osteopontin. It recognizes the sequence A-E-I-D-G-I-E-L in cytotactin.,induction: Up-regulated in small cell lung cancers (SCLC).,similarity: Belongs to the integrin alpha chain family.,similarity: Contains 7 FG-GAP repeats.,subunit: Heterodimer of an alpha and a beta subunit. Alpha-9 associates with beta-1.,tissue specificity: In airway epithelium, in the basal layer of squamous epithelium, and in smooth muscle, skeletal muscle, and hepatocytes. Abundantly expressed in fetal lung and lung cancers.,
Background	integrin subunit alpha 9(ITGA9) Homo sapiens This gene encodes an alpha integrin. Integrins are heterodimeric integral membrane glycoproteins composed of an alpha chain and a beta chain that mediate cell-cell and cell-matrix adhesion. The protein encoded by this gene, when bound to the beta 1 chain, forms an integrin that is a receptor for VCAM1, cytotactin and osteopontin. Expression of this gene has been found to be upregulated in small cell lung cancers. [provided

by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images