




STRN4 Polyclonal Antibody

Catalog No	YP-Ab-05188
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	STRN4 ZIN
Protein Name	Striatin-4 (Zinedin)
Immunogen	Synthesized peptide derived from human protein . at AA range: 250-330
Specificity	STRN4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	82kD
Cell Pathway	Cytoplasm . Membrane ; Peripheral membrane protein . Cell projection, dendritic spine . CTTNBP2-binding may regulate dendritic spine distribution. .
Tissue Specificity	Brain,Muscle,Testis,
Function	function:Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein.,miscellaneous:The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane).,similarity:Belongs to the WD repeat striatin family.,similarity:Contains 7 WD repeats.,
Background	function:Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein.,miscellaneous:The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane).,similarity:Belongs to the WD repeat striatin family.,similarity:Contains 7 WD repeats.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of Hela lysis, using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000