



# B3A2 Polyclonal Antibody

<b>Catalog No</b>	YP-Ab-05321
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse
<b>Applications</b>	WB;ELISA
<b>Gene Name</b>	SLC4A2 AE2 EPB3L1 HKB3 MPB3L
<b>Protein Name</b>	Anion exchange protein 2 (AE 2) (Anion exchanger 2) (Non-erythroid band 3-like protein) (BND3L) (Solute carrier family 4 member 2)
<b>Immunogen</b>	Synthesized peptide derived from human protein . at AA range: 830-910
<b>Specificity</b>	B3A2 Polyclonal Antibody detects endogenous levels of protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	136kD
<b>Cell Pathway</b>	Membrane; Multi-pass membrane protein.
<b>Tissue Specificity</b>	Brain,
<b>Function</b>	alternative products:Additional isoforms seem to exist,function:Plasma membrane anion exchange protein of wide distribution.,similarity:Belongs to the anion exchanger (TC 2.A.31) family.,
<b>Background</b>	alternative products:Additional isoforms seem to exist,function:Plasma membrane anion exchange protein of wide distribution.,similarity:Belongs to the anion exchanger (TC 2.A.31) family.,
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images