



CTNA2 Polyclonal Antibody

Catalog No	YP-Ab-05436
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	CTNNA2 CAPR
Protein Name	Catenin alpha-2 (Alpha N-catenin) (Alpha-catenin-related protein)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	CTNA2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	104kD
Cell Pathway	Cell membrane ; Peripheral membrane protein ; Cytoplasmic side . Cytoplasm . Cytoplasm, cytoskeleton . Cell junction, adherens junction . Cell projection, axon . Nucleus .
Tissue Specificity	Expressed in neural tissues, with strongest expression in fetal and adult brain. Expressed in the developing cortical plate and marginal zone of 20-week-old human fetal brain.
Function	similarity:Belongs to the vinculin/alpha-catenin family.,
Background	similarity:Belongs to the vinculin/alpha-catenin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images