



GCFC1 Polyclonal Antibody

Catalog No	YP-Ab-05598
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	GCFC1 C21orf66 GCFC
Protein Name	GC-rich sequence DNA-binding factor 1
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	GCFC1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	100kD
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitous.
Function	function:Possible transcription factor.,similarity:Belongs to the GCF family.,tissue specificity:Ubiquitous.,
Background	PAX3 and PAX7 binding protein 1(PAXBP1) Homo sapiens This gene encodes a protein that may bind to GC-rich DNA sequences, which suggests its involvement in the regulation of transcription. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jun 2009],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images