



LARP1 Polyclonal Antibody

Catalog No	YP-Ab-05689
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	LARP1 KIAA0731 LARP
Protein Name	La-related protein 1 (La ribonucleoprotein domain family member 1)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	LARP1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	120kD
Cell Pathway	Cytoplasm . Cytoplasmic granule . Colocalizes with RPTOR and PABPC1 in cytoplasmic granules that resemble stress granules. .
Tissue Specificity	Brain,Cervix,Cervix carcinoma,Colon carcinoma,Epithelium,Ma
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the LARP family.,similarity:Contains 1 HTH La-type RNA-binding domain.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the LARP family.,similarity:Contains 1 HTH La-type RNA-binding domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images