



# PCF11 Polyclonal Antibody

<b>Catalog No</b>	YP-Ab-05967
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Rat;Mouse;
<b>Applications</b>	WB;ELISA
<b>Gene Name</b>	PCF11 KIAA0824
<b>Protein Name</b>	Pre-mRNA cleavage complex 2 protein Pcf11 (Pre-mRNA cleavage complex II protein Pcf11)
<b>Immunogen</b>	Synthesized peptide derived from human protein . at AA range: 490-570
<b>Specificity</b>	PCF11 Polyclonal Antibody detects endogenous levels of protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	171kD
<b>Cell Pathway</b>	Nucleus.
<b>Tissue Specificity</b>	Brain,Epithelium,Testis,Uterus,
<b>Function</b>	function:Component of pre-mRNA cleavage complex II.,sequence caution:Contaminating sequence. Potential poly-A sequence.,similarity:Contains 1 CID domain.,
<b>Background</b>	The protein encoded by this gene binds to CLP1 to form pre-mRNA cleavage factor IIm. The encoded protein is necessary for efficient Pol II transcription termination and may be involved in degradation of the 3' product of polyA site cleavage. [provided by RefSeq, Oct 2016],
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images