



POTEE Polyclonal Antibody

Catalog No	YP-Ab-05978
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	POTEE A26C1A POTE2
Protein Name	POTE ankyrin domain family member E (ANKRD26-like family C member 1A) (Prostate, ovary, testis-expressed protein on chromosome 2) (POTE-2)
Immunogen	Synthesized peptide derived from human protein . at AA range: 350-430
Specificity	POTEE Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	118kD
Cell Pathway	extracellular space,extracellular exosome,blood microparticle,
Tissue Specificity	Hepatocellular carcinoma,Prostate,
Function	similarity:Belongs to the actin family.,similarity:Contains 5 ANK repeats.,similarity:In the C-terminal section; belongs to the actin family.,similarity:In the N-terminal section; belongs to the POTE family.,
Background	similarity:Belongs to the actin family.,similarity:Contains 5 ANK repeats.,similarity:In the C-terminal section; belongs to the actin family.,similarity:In the N-terminal section; belongs to the POTE family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images