



# DEN5A Polyclonal Antibody

<b>Catalog No</b>	YP-Ab-06005
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse
<b>Applications</b>	WB;ELISA
<b>Gene Name</b>	DENND5A KIAA1091 RAB6IP1
<b>Protein Name</b>	DENN domain-containing protein 5A (Rab6-interacting protein 1) (Rab6IP1)
<b>Immunogen</b>	Synthesized peptide derived from human protein . at AA range: 650-730
<b>Specificity</b>	DEN5A Polyclonal Antibody detects endogenous levels of protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	141kD
<b>Cell Pathway</b>	Golgi apparatus membrane .
<b>Tissue Specificity</b>	Brain,Placenta,Skin,Uterus,
<b>Function</b>	function:May be involved in Rab6-mediated GTPase signaling.,similarity:Belongs to the RAB6IP1 family.,similarity:Contains 1 dDENN domain.,similarity:Contains 1 DENN domain.,similarity:Contains 1 PLAT domain.,similarity:Contains 1 uDENN domain.,similarity:Contains 2 RUN domains.,subunit:Probably interacts with RAB6A bound to GTP.,
<b>Background</b>	function:May be involved in Rab6-mediated GTPase signaling.,similarity:Belongs to the RAB6IP1 family.,similarity:Contains 1 dDENN domain.,similarity:Contains 1 DENN domain.,similarity:Contains 1 PLAT domain.,similarity:Contains 1 uDENN domain.,similarity:Contains 2 RUN domains.,subunit:Probably interacts with RAB6A bound to GTP.,
<b>matters needing attention</b>	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**