



SCAFB Polyclonal Antibody

Catalog No	YP-Ab-06161
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	SCAF11 CASP11 SFRS2IP SIP1 SRSF2IP
Protein Name	Protein SCAF11 (CTD-associated SR protein 11) (Renal carcinoma antigen NY-REN-40) (SC35-interacting protein 1) (SR-related and CTD-associated factor 11) (SRSF2-interacting protein) (Serine/arginine-rich protein)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SCAFB Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	160kD
Cell Pathway	Nucleus .
Tissue Specificity	Widely expressed.
Function	function:Plays a role in pre-mRNA alternative splicing by regulating spliceosome assembly.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Contaminating sequence. Potential poly-A sequence.,subunit:Interacts with SFRS2, U2AF2 and SNRNP70.,tissue specificity:Widely expressed.,
Background	function:Plays a role in pre-mRNA alternative splicing by regulating spliceosome assembly.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Contaminating sequence. Potential poly-A sequence.,subunit:Interacts with SFRS2, U2AF2 and SNRNP70.,tissue specificity:Widely expressed.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images