



SHE Polyclonal Antibody

Catalog No	YP-Ab-06162
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	SHE
Protein Name	SH2 domain-containing adapter protein E
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SHE Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	54kD
Cell Pathway	
Tissue Specificity	Spleen,
Function	similarity:Contains 1 SH2 domain.,
Background	similarity:Contains 1 SH2 domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

