



# S6A13 Polyclonal Antibody

<b>Catalog No</b>	YP-Ab-06195
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Rat;Mouse
<b>Applications</b>	WB;ELISA
<b>Gene Name</b>	SLC6A13 GAT2
<b>Protein Name</b>	Sodium- and chloride-dependent GABA transporter 2 (GAT-2) (Solute carrier family 6 member 13)
<b>Immunogen</b>	Synthesized peptide derived from part region of human protein
<b>Specificity</b>	S6A13 Polyclonal Antibody detects endogenous levels of protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	66kD
<b>Cell Pathway</b>	Cell membrane ; Multi-pass membrane protein .
<b>Tissue Specificity</b>	Expressed in brain, kidney, lung, liver and testis.
<b>Function</b>	function:Terminates the action of GABA by its high affinity sodium-dependent reuptake into presynaptic terminals.,similarity:Belongs to the sodium:neurotransmitter symporter (SNF) family.,tissue specificity:Expressed in brain, kidney, lung, liver and testis.,
<b>Background</b>	function:Terminates the action of GABA by its high affinity sodium-dependent reuptake into presynaptic terminals.,similarity:Belongs to the sodium:neurotransmitter symporter (SNF) family.,tissue specificity:Expressed in brain, kidney, lung, liver and testis.,
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images