



NAC2 Polyclonal Antibody

Catalog No	YP-Ab-06204
Isotype	IgG
Reactivity	Human;Rat
Applications	WB;ELISA
Gene Name	SLC8A2 KIAA1087 NCX2
Protein Name	Sodium/calcium exchanger 2 (Na ⁺)/Ca ²⁺ -exchange protein 2)
Immunogen	Synthesized peptide derived from human protein . at AA range: 380-460
Specificity	NAC2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	101kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Basolateral cell membrane ; Multi-pass membrane protein . Perikaryon . Cell projection, dendrite . Cell projection, dendritic spine .
Tissue Specificity	Amygdala,Brain,
Function	enzyme regulation:By ATP.,function:Rapidly transports Ca ²⁺ during excitation-contraction coupling. Ca ²⁺ is extruded from the cell during relaxation so as to prevent overloading of intracellular stores.,similarity:Belongs to the sodium/potassium/calcium exchanger family. SLC8 subfamily.,similarity:Contains 2 Calx-beta domains.,
Background	enzyme regulation:By ATP.,function:Rapidly transports Ca ²⁺ during excitation-contraction coupling. Ca ²⁺ is extruded from the cell during relaxation so as to prevent overloading of intracellular stores.,similarity:Belongs to the sodium/potassium/calcium exchanger family. SLC8 subfamily.,similarity:Contains 2 Calx-beta domains.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images