



SL9C2 Polyclonal Antibody

Catalog No	YP-Ab-06208
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	SLC9C2 SLC9A11
Protein Name	Sodium/hydrogen exchanger 11 (Na ⁺)/H ⁺ exchanger 11) (NHE-11) (Solute carrier family 9 member 11) (Solute carrier family 9 member C2)
Immunogen	Synthesized peptide derived from part region of human protein. AA range 10-60
Specificity	SL9C2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	123kD
Cell Pathway	Membrane; Multi-pass membrane protein.
Tissue Specificity	Testis,
Function	domain:Contains an ion transport-like region is related to the membrane segments of voltage-gated ion channels.,function:Involved in pH regulation.,similarity:Belongs to the monovalent cation:proton antiporter 1 (CPA1) transporter (TC 2.A.36) family.,similarity:Contains 1 cyclic nucleotide-binding domain.,
Background	domain:Contains an ion transport-like region is related to the membrane segments of voltage-gated ion channels.,function:Involved in pH regulation.,similarity:Belongs to the monovalent cation:proton antiporter 1 (CPA1) transporter (TC 2.A.36) family.,similarity:Contains 1 cyclic nucleotide-binding domain.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images