



SPEG Polyclonal Antibody

Catalog No	YP-Ab-06241
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	IHC;IF
Gene Name	SPEG APEG1 KIAA1297
Protein Name	Striated muscle preferentially expressed protein kinase (EC 2.7.11.1) (Aortic preferentially expressed protein 1) (APEG-1)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SPEG Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC-p 1:50-300. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	359kD
Cell Pathway	[Isoform 3]: Nucleus.
Tissue Specificity	Isoform 1 is preferentially expressed in striated muscle. Non-kinase form such as isoform 3 is predominantly expressed in the aorta. Isoform 3 appears to be expressed only in highly differentiated ASMC in normal vessel walls and down-regulated in dedifferentiated ASMC in vivo. In response to vascular injuries ASMC dedifferentiate and change from a quiescent and contractile phenotype to a proliferative and synthetic phenotype. This proliferation of vascular smooth muscle cells is one of the most prominent features of atherosclerosis.
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,function:Isoform 3 may have a role in regulating the growth and differentiation of arterial smooth muscle cells.,induction:Isoform 3 is quickly down-regulated in response to vascular injury, when ASMC cells change from a quiescent to a proliferative phenotype.,miscellaneous:Expression is under the tight control of the locus control region (LCRs).,PTM:May be autophosphorylated. Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family.,similarity:Contains 2 fibronectin type-III domains.,similarity:Contains 2 protein kinase domains.,similarity:Contains 9 Ig-like (immunoglobulin-like) domains.,subunit:Isoform 3 is found as a monomer



or homodimer.,tissue specificity:Isoform 1 is preferentially expressed in striated muscle. Non-kinase form such as

Background

This gene encodes a protein with similarity to members of the myosin light chain kinase family. This protein family is required for myocyte cytoskeletal development. Along with the desmin gene, expression of this gene may be controlled by the desmin locus control region. Mutations in this gene are associated with centronuclear myopathy 5. [provided by RefSeq, Jun 2016],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Tris-EDTA,pH9.0 was used for antigen retrieval. 2 Antibody was diluted at 1:200(4° overnight.3,Secondary antibody was diluted at 1:200(room temperature, 45min).