



UBP29 Polyclonal Antibody

Catalog No	YP-Ab-06339
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	USP29
Protein Name	Ubiquitin carboxyl-terminal hydrolase 29 (EC 3.4.19.12) (Deubiquitinating enzyme 29) (Ubiquitin thioesterase 29) (Ubiquitin-specific-processing protease 29)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	UBP29 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	101kD
Cell Pathway	Cytoplasm, perinuclear region . Localizes to perinuclear region in response to herpes simplex virus-1 (HSV-1) infection. .
Tissue Specificity	
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol.,similarity:Belongs to the peptidase C19 family.,
Background	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol.,similarity:Belongs to the peptidase C19 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

