



UBP3 Polyclonal Antibody

Catalog No	YP-Ab-06340
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	USP3
Protein Name	Ubiquitin carboxyl-terminal hydrolase 3 (EC 3.4.19.12) (Deubiquitinating enzyme 3) (Ubiquitin thioesterase 3) (Ubiquitin-specific-processing protease 3)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	UBP3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	57kD
Cell Pathway	Nucleus . Localizes preferentially with monoubiquitinated H2A to chromatin.
Tissue Specificity	Expressed in all tissues examined, with strongest expression in pancreas.
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol.,similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 UBQ-type zinc finger.,tissue specificity:Expressed in all tissues examined, with strongest expression in pancreas.,
Background	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol.,similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 UBQ-type zinc finger.,tissue specificity:Expressed in all tissues examined, with strongest expression in pancreas.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images