




ARAID Polyclonal Antibody

Catalog No	YP-Ab-06423
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	ATRAID APR3 C2orf28 HSPC013 UNQ214/PRO240
Protein Name	All-trans retinoic acid-induced differentiation factor (Apoptosis-related protein 3) (APR-3) (p18)
Immunogen	Synthesized peptide derived from human protein . at AA range: 1-80
Specificity	ARAID Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	25kD
Cell Pathway	Nucleus envelope . Cell membrane ; Single-pass membrane protein . Colocalizes with NELL1 on the nuclear envelope and the perinuclear region (PubMed:21723284).
Tissue Specificity	Weakly expressed in hematopoietic cell lines.
Function	sequence caution:Translation N-terminally shortened.,similarity:Contains 1 EGF-like domain.,tissue specificity:Weakly expressed in hematopoietic cell lines.,
Background	This gene is thought to be involved in apoptosis, and may also be involved in hematopoietic development and differentiation. The use of alternative splice sites and promoters result in multiple transcript variants encoding different isoforms.[provided by RefSeq, Dec 2009],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Jarkat cells,
primary antibody was diluted at 1:1000, 4°over night