



TAF1B Polyclonal Antibody

Catalog No	YP-Ab-06498
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	TAF1B
Protein Name	TATA box-binding protein-associated factor RNA polymerase I subunit B (RNA polymerase I-specific TBP-associated factor 63 kDa) (TAFI63) (TATA box-binding protein-associated factor 1B) (TBP-associated
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TAF1B Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	64kD
Cell Pathway	Nucleus, nucleolus.
Tissue Specificity	Colon,Corpus callosum,
Function	function:Component of the transcription factor SL1/TIFIB complex involved in the assembly of the PIC (preinitiation complex).,similarity:Belongs to the TAF1 family.,subunit:Component of the transcription factor SL1/TIFIB complex which is composed of TBP and at least TAF1A, TAF1B and TAF1C. In the complex, interacts directly with TBP, TAF1A and TAF1C. Interaction of the SL1/TIFIB subunits with TBP excludes interaction with the transcription factor IID (TFIID) subunits. This mutually exclusive binding directs the formation of promoter and RNA polymerase-selective TBP-TAF complexes.,
Background	Initiation of transcription by RNA polymerase I requires the formation of a complex composed of the TATA-binding protein (TBP) and three TBP-associated factors (TAFs) specific for RNA polymerase I. This complex, known as SL1, binds to the core promoter of ribosomal RNA genes to position the polymerase properly and acts as a channel for regulatory signals. This gene encodes one of the SL1-specific TAFs. Multiple transcript variants encoding different isoforms have

been found for this gene. [provided by RefSeq, Jan 2016],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images