



B2L13 Polyclonal Antibody

Catalog No	YP-Ab-06722
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	BCL2L13 MIL1 CD003
Protein Name	Bcl-2-like protein 13 (Bcl2-L-13) (Bcl-rambo) (Protein Mil1)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	B2L13 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	53kD
Cell Pathway	[Isoform 2]: Mitochondrion membrane ; Single-pass membrane protein . Nucleus .; [Isoform 1]: Nucleus .
Tissue Specificity	Ubiquitous, with the highest levels of expression in heart, placenta and pancreas.
Function	alternative products:Experimental confirmation may be lacking for some isoforms,function:May promote the activation of caspase-3 and apoptosis.,similarity:Belongs to the Bcl-2 family.,subunit:Monomer.,tissue specificity:Ubiquitous, with the highest levels of expression in heart, placenta and pancreas.,
Background	This gene encodes a mitochondrially-localized protein with conserved B-cell lymphoma 2 homology motifs. Overexpression of the encoded protein results in apoptosis. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jul 2012],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images