




ACTBM Polyclonal Antibody

Catalog No	YP-Ab-06775
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	POTEKP ACTBL3 FKSG30
Protein Name	Putative beta-actin-like protein 3 (Kappa-actin) (POTE ankyrin domain family member K)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ACTBM Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	41kD
Cell Pathway	Cytoplasm, cytoskeleton .
Tissue Specificity	Expressed in some hepatocellular carcinomas.
Function	similarity:Belongs to the actin family.,subunit:Interacts with PFN1 and PFDN1. Does not interact with PFN2.,tissue specificity:Expressed in some hepatocellular carcinomas.,
Background	similarity:Belongs to the actin family.,subunit:Interacts with PFN1 and PFDN1. Does not interact with PFN2.,tissue specificity:Expressed in some hepatocellular carcinomas.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4°over night