



LRIG3 Polyclonal Antibody

Catalog No	YP-Ab-06836
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	LRIG3 LIG3 UNQ287/PRO326/PRO335
Protein Name	Leucine-rich repeats and immunoglobulin-like domains protein 3 (LIG-3)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	LRIG3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	123kD
Cell Pathway	Cell membrane ; Single-pass type I membrane protein . Cytoplasmic vesicle membrane ; Single-pass type I membrane protein . Detected in cytoplasmic vesicles when coexpressed with ERBB4. .
Tissue Specificity	Widely expressed.
Function	similarity:Contains 17 LRR (leucine-rich) repeats.,similarity:Contains 3 Ig-like C2-type (immunoglobulin-like) domains.,tissue specificity:Widely expressed.,
Background	similarity:Contains 17 LRR (leucine-rich) repeats.,similarity:Contains 3 Ig-like C2-type (immunoglobulin-like) domains.,tissue specificity:Widely expressed.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

