



AFAM Polyclonal Antibody

Catalog No	YP-Ab-06939
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	AFM ALB2 ALBA
Protein Name	Afamin (Alpha-albumin) (Alpha-Alb)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	AFAM Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	65kD
Cell Pathway	Secreted .
Tissue Specificity	High level detected in plasma but also in extravascular fluids such as follicular and cerebrospinal fluids (at protein level).
Function	function:Possible role in the transport of yet unknown ligand.,PTM:N-glycosylated.,similarity:Belongs to the ALB/AFP/VDB family.,similarity:Contains 3 albumin domains.,
Background	This gene is a member of the albumin gene family, which is comprised of four genes that localize to chromosome 4 in a tandem arrangement. These four genes encode structurally-related serum transport proteins that are known to be evolutionarily related. The protein encoded by this gene is regulated developmentally, expressed in the liver and secreted into the bloodstream. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images