




DPPA4 Polyclonal Antibody

Catalog No	YP-Ab-06944
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	DPPA4
Protein Name	Developmental pluripotency-associated protein 4
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	DPPA4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	33kD
Cell Pathway	Nucleus . Associated with transcriptionally active chromatin. .
Tissue Specificity	Teratocarcinoma,Testis,
Function	function:May play a role in maintaining cell pluripotentiality.,
Background	function:May play a role in maintaining cell pluripotentiality.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night