




KLD8B rabbit pAb

Catalog No	YP-Ab-07942
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	KLHDC8B FP17659
Protein Name	KLD8B
Immunogen	Synthesized peptide derived from human KLD8B AA range: 120-170
Specificity	This antibody detects endogenous levels of KLD8B at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.56% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Kelch domain-containing protein 8B
Observed Band	38kD
Cell Pathway	Cytoplasm . Midbody . In mitotic cells, concentrates in the midbody of the cytoplasmic bridge linking daughter cells as they are about to separate during cytokinesis. .
Tissue Specificity	Adrenal cortex, Brain,
Function	similarity: Contains 8 Kelch repeats.,
Background	This gene encodes a protein which forms a distinct beta-propeller protein structure of kelch domains allowing for protein-protein interactions. Mutations in this gene have been associated with Hodgkin lymphoma. [provided by RefSeq, Sep 2010],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC-12 cells,
primary antibody was diluted at 1:1000, 4°over night