




F162A rabbit pAb

Catalog No	YP-Ab-08364
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	FAM162A C3orf28 E2IG5 DC16 FWP001
Protein Name	F162A
Immunogen	Synthesized peptide derived from human F162A AA range: 78-128
Specificity	This antibody detects endogenous levels of F162A at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Mitochondrion membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	induction:By 17-beta-estradiol.,similarity:Belongs to the UPF0389 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night