




GABP2 rabbit pAb

Catalog No	YP-Ab-08434
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	GABPB2
Protein Name	GABP2
Immunogen	Synthesized peptide derived from human GABP2 AA range: 314-364
Specificity	This antibody detects endogenous levels of GABP2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:May function as transcription factor capable of interacting with purine rich repeats (GA repeats).,similarity:Contains 5 ANK repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night