




PUS3 rabbit pAb

Catalog No	YP-Ab-08532
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PUS3 FKSG32
Protein Name	PUS3
Immunogen	Synthesized peptide derived from human PUS3 AA range: 335-385
Specificity	This antibody detects endogenous levels of PUS3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	catalytic activity:tRNA uridine = tRNA pseudouridine.,function:Formation of pseudouridine at position 39 in the anticodon stem and loop of transfer RNAs.,similarity:Belongs to the tRNA pseudouridine synthase truA family.,
Background	The protein encoded by this gene catalyzes the formation of tRNA pseudouridine from tRNA uridine at position 39 in the anticodon stem and loop of transfer RNAs. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2012],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night