




LRC69 rabbit pAb

Catalog No	YP-Ab-08604
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	LRRC69
Protein Name	LRC69
Immunogen	Synthesized peptide derived from human LRC69 AA range: 142-192
Specificity	This antibody detects endogenous levels of LRC69 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the LRRC69 family.,similarity:Contains 9 LRR (leucine-rich) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from DU145 cells,
primary antibody was diluted at 1:1000, 4° over night