




MSD1 rabbit pAb

Catalog No	YP-Ab-08679
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	MSANTD1 C4orf44
Protein Name	MSD1
Immunogen	Synthesized peptide derived from human MSD1 AA range: 61-111
Specificity	This antibody detects endogenous levels of MSD1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	nuclear body,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night