




CK096 rabbit pAb

Catalog No	YP-Ab-08734
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	C11orf96 AG2
Protein Name	CK096
Immunogen	Synthesized peptide derived from human CK096. AA rang 1-50
Specificity	This antibody detects endogenous levels of CK096 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	PTM: Phosphorylated upon DNA damage, probably by ATM or ATR.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night