




OSBP2 rabbit pAb

Catalog No	YP-Ab-08736
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	OSBP2 KIAA1664 ORP4 OSBPL4
Protein Name	OSBP2
Immunogen	Synthesized peptide derived from human OSBP2 AA range: 86-136
Specificity	This antibody detects endogenous levels of OSBP2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane; Peripheral membrane protein.
Tissue Specificity	Expressed mainly in retina, testis, and fetal liver.
Function	function: Binds 7-ketocholesterol., similarity: Belongs to the OSBP family., similarity: Contains 1 PH domain., tissue specificity: Expressed mainly in retina, testis, and fetal liver.,
Background	The protein encoded by this gene contains a pleckstrin homology (PH) domain and an oxysterol-binding region. It binds oxysterols such as 7-ketocholesterol and may inhibit their cytotoxicity. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Sep 2013],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4°over night