




AMER3 rabbit pAb

Catalog No	YP-Ab-08750
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM123C
Protein Name	AMER3
Immunogen	Synthesized peptide derived from human AMER3 AA range: 676-726
Specificity	This antibody detects endogenous levels of AMER3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane ; Peripheral membrane protein . Translocates to the cell membrane following binding to PtdIns(4,5)P2. .
Tissue Specificity	
Function	similarity:Belongs to the FAM123 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SH-SY5Y cells,
primary antibody was diluted at 1:1000, 4° over night