




# ZDH12 rabbit pAb

|                                  |                                                                                                                                                                                                                                                                         |
|----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                | YP-Ab-08775                                                                                                                                                                                                                                                             |
| <b>Isotype</b>                   | IgG                                                                                                                                                                                                                                                                     |
| <b>Reactivity</b>                | Human; Mouse; Rat                                                                                                                                                                                                                                                       |
| <b>Applications</b>              | WB                                                                                                                                                                                                                                                                      |
| <b>Gene Name</b>                 | ZDHHC12 ZNF400 PSEC0008                                                                                                                                                                                                                                                 |
| <b>Protein Name</b>              | ZDH12                                                                                                                                                                                                                                                                   |
| <b>Immunogen</b>                 | Synthesized peptide derived from human ZDH12 AA range: 83-133                                                                                                                                                                                                           |
| <b>Specificity</b>               | This antibody detects endogenous levels of ZDH12 at Human/Mouse/Rat                                                                                                                                                                                                     |
| <b>Formulation</b>               | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                 |
| <b>Source</b>                    | Polyclonal, Rabbit, IgG                                                                                                                                                                                                                                                 |
| <b>Purification</b>              | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                                                                                                                               |
| <b>Dilution</b>                  | WB 1: 500-2000                                                                                                                                                                                                                                                          |
| <b>Concentration</b>             | 1 mg/ml                                                                                                                                                                                                                                                                 |
| <b>Purity</b>                    | ≥90%                                                                                                                                                                                                                                                                    |
| <b>Storage Stability</b>         | -20°C/1 year                                                                                                                                                                                                                                                            |
| <b>Synonyms</b>                  |                                                                                                                                                                                                                                                                         |
| <b>Observed Band</b>             |                                                                                                                                                                                                                                                                         |
| <b>Cell Pathway</b>              | Golgi apparatus membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein .                                                                                                                                                 |
| <b>Tissue Specificity</b>        | Widely expressed.                                                                                                                                                                                                                                                       |
| <b>Function</b>                  | catalytic activity: Palmitoyl-CoA + protein-cysteine = S-palmitoyl protein + CoA.; domain: The DHHC domain is required for palmitoyltransferase activity.; similarity: Belongs to the DHHC palmitoyltransferase family.; similarity: Contains 1 DHHC-type zinc finger.; |
| <b>Background</b>                |                                                                                                                                                                                                                                                                         |
| <b>matters needing attention</b> | Avoid repeated freezing and thawing!                                                                                                                                                                                                                                    |
| <b>Usage suggestions</b>         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                                                                                                                                       |



## Products Images



Western blot analysis of lysates from A431 cells,  
primary antibody was diluted at 1:1000, 4°over night