




TSACC rabbit pAb

Catalog No	YP-Ab-08836
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TSACC C1orf182
Protein Name	TSACC
Immunogen	Synthesized peptide derived from human TSACC AA range: 57-107
Specificity	This antibody detects endogenous levels of TSACC at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	cytoplasm,
Tissue Specificity	Expressed in testis but is absent from mature sperm.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night