




SRPX rabbit pAb

Catalog No	YP-Ab-08937
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	SRPX ETX1
Protein Name	SRPX
Immunogen	Synthesized peptide derived from human SRPX AA range: 234-284
Specificity	This antibody detects endogenous levels of SRPX at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1: 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell surface . Possibly surface of photoreceptor cell.
Tissue Specificity	Retina and heart; less in placenta, pancreas, lung, liver, skeletal muscle, kidney and brain.
Function	alternative products: Additional isoforms seem to exist, function: May be involved in phagocytosis during disk shedding, cell adhesion to cells other than the pigment epithelium or signal transduction., similarity: Contains 1 HYR domain., similarity: Contains 3 Sushi (CCP/SCR) domains., subcellular location: Possibly surface of photoreceptor cell., tissue specificity: Retina and heart; less in placenta, pancreas, lung, liver, skeletal muscle, kidney and brain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night