



K1407 rabbit pAb

Catalog No	YP-Ab-09186
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	KIAA1407
Protein Name	Coiled-coil domain-containing protein KIAA1407
Immunogen	Synthesized peptide derived from human K1407 AA range: 720-770
Specificity	This antibody detects endogenous levels of human K1407
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

