



YN4361

UN93A rabbit pAb

Catalog No	YP-Ab-09201
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	UNC93A
Protein Name	Protein unc-93 homolog A (HmUnc-93A) (Unc-93A)
Immunogen	Synthesized peptide derived from human UN93A AA range: 312-362
Specificity	This antibody detects endogenous levels of human UN93A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in testis, small intestine, spleen, prostate and ovary.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images