



LOXE3 rabbit pAb

Catalog No	YP-Ab-09203
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ALOXE3
Protein Name	Epidermis-type lipxygenase 3 (e-LOX-3) (EC 1.13.11.-)
Immunogen	Synthesized peptide derived from human LOXE3 AA range: 467-517
Specificity	This antibody detects endogenous levels of human LOXE3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
molecular weight	78kD
Observed Band	
Cell Pathway	Cytoplasm .
Tissue Specificity	Predominantly expressed in skin.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

